

S100A14-S100 calcium binding protein A14 Gene

Reference: PTXBC005019

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BCMP84; S100A15; protein S100-A14; S114; breast cancer membrane protein 84; S100 calcium binding protein A14		

For research use only.

GENE DESCRIPTION	S100 calcium binding protein A14
GENE ID	57402
ORIGIN SPECIES	Human
ACCESSIONS	BC005019
VECTOR	pDONR223
NCBI SYMBOL	S100A14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacagtgtcgtcagccaacgcagaggatgctcaggaattcagtgatgtggagaggccattgagaccctcatcaagaacttcaccagtactcgtggagggtgggaaggagacgtgaccccttctgagctacgggacctggtcaccagcagctgccccatctatgccgag
caactgtggcctggaagagaaaattgccaacctgggcagctgcaatgactctaaactggagttcaggagtttctgggagctgattggagaagcggccaagagtgtagagctggagaggcctgtccggggcactga