

FUSIP1-FUS interacting protein (serine/arginine-rich) 1 Gene

Reference: PTXBC005039

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FUSIP1; FUSIP2; NSSR; PPP1R149; SFRS13; SFRS13A; SRp38; SRp40; TASR; TASR1; TASR2; serine/arginine-rich splicing factor 10; 40 kDa SR-repressor protein; FUS interacting protein (serine-arginine rich) 1; FUS-interacting protein (serine-arginine rich) 2; S		

For research use only.

GENE DESCRIPTION	FUS interacting protein (serine/arginine-rich) 1
GENE ID	10772
ORIGIN SPECIES	Human
ACCESSIONS	BC005039
VECTOR	pDONR223
NCBI SYMBOL	FUSIP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcccgctacctgctcccccaacacgtctctgttcaggaacgtggccgacgacaccaggtctgaagactgcggcgtgaatttggtcgttatggtcctatagttgatgtgtatgttccacttgatttctacactcgcgcgccaagagattgcttatgtcaatttgaggatgttcgtgat
gctgaagacgtttacataatttggacagaaagtggatttggacggcagattgaaatcacgtttgccagggggatcgaagacaccaaatacagatgaaagccaaggaaggaggaatgtgtacagttcttcacgtatgatgattatgacagatcacagcgttctagaagccgaagtat
gaaaggaggagatcaagaagtcggcttttattacaactatagaagatcgtatagtcctagaacagtagaccgactggaagaccacgcgtgacgaagccattccgacaatgatagattcaaacaccgaaatcgatcttttcaagatctaaatccaattcaagatcacggtccaagtcc
cagcccaagaaagaatgaaggctaaatcacgtttaggtctgcatctcaccaaaactagagcaccctctaaaacagattccaaacacattataagtcgtcgaagatgaaaggaatcaaggaaaaagaaccacctagatccaaatctcagtcagatcacagtcaggtctag
gtcaaaatctagatcaaggctcttgactagtcctaaagtccagtgccactga