

MGC31957-hypothetical protein MGC31957 Gene

Reference: PTXBC005043

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC31957
GENE ID	254896
ORIGIN SPECIES	Human
ACCESSIONS	BC005043
VECTOR	pDONR223
NCBI SYMBOL	MGC31957
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaagggtgaaggagcgttctaccgttagggaaactctggggacagagcgcggcgccgctgatggccgagggcaggggtgcgaccaggaccaggacggcgctgggaaccataccatggcccgatcccaagaccctaaagttcgtcgtcatcgtcgggtcctgct
gccagtgcgtcccgccgctcctggctggggaaagagcgcacctggcgccgggagggggcaggagacggggacacggcagggatgcctggccctgggtcacctgcggccggcatgtccgggcaggacgaactcgcgtcggagtcagggagaactgggtcccggtg
ctgggcaggagggacccggccgagggagcagagaggggtccctcgtgctgccccgagcccggaaggagggaagttccagaatcgagaggggagggagtcagggtggaacccatag