



Reference: PTXBC005046

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	V75; VNT; S-protein; complement S-protein; epibolin; serum spreading factor; somatomedin B		

For research use only.

GENE DESCRIPTION	vitronectin
GENE ID	7448
ORIGIN SPECIES	Human
ACCESSIONS	BC005046
VECTOR	pDONR223
NCBI SYMBOL	VTN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr