

FOXP1-forkhead box P1 Gene

Reference: PTXBC005055

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	12CC4; HSPC215; MFH; QRF1; hFKH1B; forkhead box protein P1; fork head-related protein like B; glutamine-rich factor 1; mac-1-regulated forkhead; forkhead box P1		

For research use only.

GENE DESCRIPTION	forkhead box P1
GENE ID	27086
ORIGIN SPECIES	Human
ACCESSIONS	BC005055
VECTOR	pDONR223
NCBI SYMBOL	FOXP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgcaagaatctgggactgagacaaaaagtaacgggtcagccatccagaatgggtcggcgccagcaaccacttagagtgccggcgtctctgggagggcggtccaacggagagacgccggccgtggacatcggggcagctgacctgcccacgcccagcagcag
caacagtgccatcataaacatcagccctctaggagtcaccagcagttgcttaagagactaattcaagcccttgggagttggaagtcctgcaggtcccttgtgggagcagttgctgagacgaagatgagtgacctgtgtgtcagcctaacccttccccatttga