

HIATL2-hippocampus abundant transcript-like 2 Gene

Reference: PTXBC005058

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hippocampus abundant transcript-like 2
GENE ID	84278
ORIGIN SPECIES	Human
ACCESSIONS	BC005058
VECTOR	pDONR223
NCBI SYMBOL	HIATL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgactgttctacatgaacattttctcaacacattctcatgaatggtctcatcaagggtgaaagggcctgctctctttttgagtgccccactcattggtgccctgtctgatgtgtgggggaggaagccctttctcctcggaactgtattctttacctgttcccaatcccactgatgaggatc
agcccatgccgggtgtgtggtgcctgtacacctgtatgccagctacctgtggcaggagaatggcttga