

FKSG24-hypothetical protein MGC12972 Gene

Reference: PTXBC005064

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FKSG24; mpv17-like protein 2; MPV17 mitochondrial membrane protein-like 2; MPV17 mitochondrial inner membrane protein like 2		

For research use only.

GENE DESCRIPTION	hypothetical protein MGC12972
GENE ID	84769
ORIGIN SPECIES	Human
ACCESSIONS	BC005064
VECTOR	pDONR223
NCBI SYMBOL	FKSG24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgggtggctggcgccgctacgcgcctgtatccgccccgcagctctatccaggccgcgcgtctctcactaacacgctggcgccgctcatggcgccggtgatggcgctgcgcagctcctgggagatccgcgccggccccgccaggttttcgacca
cggcgctccgcgagcatgtttgcgggtggcgctgagcatgggtccctctcgcactactggtactgtcgtgacccgctattccctgcgtcggcctccgaggcttcccaaatgtcctcaagaaggtcctcgtggatcagctggtagcctctccattgctggcgctctgggaattctacaagg
cagactggtgcgtgtggcctgctgcgcagttcgtgaacttctctcgtgccccccaatttcgagtcacctacatcaacggcctgacgctggcgctgggacacgtacgtctcacttgaagtacggagcccagttcctctgacacccccaggctgtgtggcctggacacccgagcagac
tga