

ERC1-ELKS/RAB6-interacting/CAST family member 1 Gene

Reference: PTXBC005065

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Cast2; ELKS; ERC-1; RAB6IP2; ELKS/Rab6-interacting/CAST family member 1; RAB6 interacting protein 2		

For research use only.

GENE DESCRIPTION	ELKS/RAB6-interacting/CAST family member 1
GENE ID	23085
ORIGIN SPECIES	Human
ACCESSIONS	BC005065
VECTOR	pDONR223
NCBI SYMBOL	ERC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctaaacaagtggctggcacggcactggccagctgctctccaagcaggtggccagatcccccacgtggactttctcaggtgcagcgcctgccactctcagccactgggtgtgtgactctctcttcctcagcattctccacttccctccagaaaaacggatgga
agggaagccctctgtgacactgttctgagaagaagcatttcgggacccgatacatctgtctgtgtctgtgaacagcaaggaatctcttga