

CCDC93-coiled-coil domain containing 93 Gene

Reference: PTXBC005078

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 93; coiled-coil domain containing 93		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 93
GENE ID	54520
ORIGIN SPECIES	Human
ACCESSIONS	BC005078
VECTOR	pDONR223
NCBI SYMBOL	CCDC93
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccacgccagccttacatccgtcccatggcctgcccagtgccctgcgcaacaaaagtaaggcctgtcttaccacactctcccagcatctag