



Reference: PTXBC005090

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	agmatinase, mitochondrial; AUH; agmatine ureohydrolase (agmatinase)		
For research use only.			

GENE DESCRIPTION	agmatine ureohydrolase (agmatinase)
GENE ID	79814
ORIGIN SPECIES	Human
ACCESSIONS	BC005090
VECTOR	pDONR223
NCBI SYMBOL	AGMAT
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1