



Reference: PTXBC005110

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RGS-PX2; SCAR20; sorting nexin-14; sorting nexin 14		
For research use only.			

GENE DESCRIPTION	sorting nexin 14
GENE ID	57231
ORIGIN SPECIES	Human
ACCESSIONS	BC005110
VECTOR	pDONR223
NCBI SYMBOL	SNX14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]