

## ADFP-adipose differentiation-related protein Gene

Reference: PTXBC005127

|                |                                                                                            |                     |                |
|----------------|--------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                        | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | ADFP; ADRP; perilipin-2; adipophilin; adipose differentiation-related protein; perilipin 2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | adipose differentiation-related protein                                                          |
| <b>GENE ID</b>          | 123                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC005127                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | ADFP                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcatcctgtgcagttgatccacaaccgagtggtgactcgggtggtcaacctgcccctgtgtgagctccacgtatgacctcatgtcctcagcctatctcagtaaaaggaccagatccctacctgaagtctgtgtgagatggcagagaacgggtggaagaccatcacctccgtggcc  
atgaccagtgctctgcccacatccagaagctagagccgcaaatgacgttgccaatcacctatgctgtgaagggctagacaggattgaggagagactgctattctgaatcagccatcaactcagattgttgccaatgccaaggcgtgtgactggggcaaaagatgctgtgacgactac  
tgtgactggggccaaggtctgtgcccagcagatcacaggggtgatggacaagaccaaaaggggagtgactggcagtggtgagaagaccaagtcgtgtcagtgccagcattaacacagctctggggagtcggatgatgcagctcgtgagcagtggtgtagaaatgcactcacc  
aaatcagagctgtgtgtagaacagtagctccctcactgagggaagaaactagaagaaagcaaaaaagtggaaggattgactctggttcagaaagccaagttattatgttagactgggacccctgtctaccaagcttcactccctgctcctaccagcaggctctcagcagggttaagaagcta  
agcaaaaaagccaacagaccatttctcagctccattctactgttcacctgattgaattgccaggaagaatgtgtatagtgccaatcagaaaaatcaggatgctcaggataagctctacctctcagggtagagtggaagggagcattggatgatgatactgatgagtcacctgtgtgag  
cacattgagtcagctactctgcaattgcccgaacctgactcagcagctccagaccagtgccacacctcctgtccaacatccaaggtgtaccacagaacatccaagatcaagccaagcacatgggggtgatggcaggcagacatcactcagtggtccgcaatgctcctcttaaga  
agtgctgacagcctcctcactctagcaaggggagcagctgcagaaaatgaaggaatctttatgacgtgatgattatctgttaacaacacgcccctcaactggctgtaggtccctttatcctcagctgactgagctcagaatgctcaggaccaagggtgcagagatggacaagagcag  
ccaggagaccagcagctgagcataaaactcattaa