

SLC35A3-solute carrier family 35 (UDP-N-acetylglucosamine (UDP-GlcNAc) transporter), member A3 Gene

Reference: PTXBC005136

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AMRS; UDP-N-acetylglucosamine transporter; golgi UDP-GlcNAc transporter; solute carrier family 35 (UDP-N-acetylglucosamine (UDP-GlcNAc) transporter), member 3; solute carrier family 35 (UDP-N-acetylglucosamine (UDP-GlcNAc) transporter), member A3; solute		

For research use only.

GENE DESCRIPTION	solute carrier family 35 (UDP-N-acetylglucosamine (UDP-GlcNAc) transporter), member A3
GENE ID	23443
ORIGIN SPECIES	Human
ACCESSIONS	BC005136
VECTOR	pDONR223
NCBI SYMBOL	SLC35A3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcgccaacctaaaatacgtttccctgggaatttggcttccagactaccagtttgggttaacaatgcgttatccagaactttaaagaagaagacccgttatctatctctacagcagtggttggctgaactttgaagataatggcctgcatatttgggtctacaagacagcaaatgta
gtctaaagacactgaatcgagtactacatgatgaaattcttaataaacctatggaacacctaacttgcctatccatcagggatctatctctcagaataattactgtatgtggcactatcaaatctagatgcagctactatcagggtacgtatcagttgaaaattctacaacagcattatctctgt
gtctatgcttagtaaaaaattgggtgtataccagtggtgctccctagtaatttgatgacaggagttgctttgtacagtgccctcagattcagcttgattctaaggaaacttcagctgggtctcaattttaggactcatggcagttctcacagcatgttttcaagtggttgggtgttactttg
agaaaattctaaagaacaaacaatcagtggtgataagaatattcagctgtgtcttttcttgagccatcctttaa