

CALM3-calmodulin 3 (phosphorylase kinase, delta) Gene

Reference: PTXBC005137

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CaM; CaMIII; HEL-S-72; PHKD; PHKD3; calmodulin; epididymis secretory protein Li 72; phosphorylase kinase subunit delta; prepro-calmodulin 3; calmodulin 3		

For research use only.

GENE DESCRIPTION	calmodulin 3 (phosphorylase kinase, delta)
GENE ID	808
ORIGIN SPECIES	Human
ACCESSIONS	BC005137
VECTOR	pDONR223
NCBI SYMBOL	CALM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaccagctgactgaggagcagattgcagagttcaaggaggccttctcctctttgacaaggatggagatggcactatcaccaccaaggagttgggacagtgatgagatccctgggacagaacccactgaagcagagctgcaggatgatcaatgagggtgatgcagatgg
gaacgggaccattgacttcccgagttcctgacctgatggccagaagatgaaggacacagacagtgaggaggagatccgagagcggttcctgtctttgacaaggatgggaatggctacatcagcgccgcagagctgcgtcacgtaatgacgaacctggggagaagctgaccga
tgaggaggtggatgagatgatcaggaggctgacatcgatggagatggccaggtcaattatgaagagttgtacagatgatgactgcaaagtga