

ARL6IP5-ADP-ribosylation-like factor 6 interacting protein 5 Gene

Reference: PTXBC005143

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DERP11; GTRAP3-18; HSPC127; JWA; PRAF3; Yip6b; addicisin; hp22; jmx; PRA1 family protein 3; ADP-ribosylation factor GTPase 6 interacting protein 5; ADP-ribosylation factor-like 6 interacting protein 5; ADP-ribosylation factor-like protein 6-interacting pro		

For research use only.

GENE DESCRIPTION	ADP-ribosylation-like factor 6 interacting protein 5
GENE ID	10550
ORIGIN SPECIES	Human
ACCESSIONS	BC005143
VECTOR	pDONR223
NCBI SYMBOL	ARL6IP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgttaatatcgcctccgctgggacgattctcccggttccgacgtttgcccggccgacgttcaggacattccaaatggaacaaccgcgtagtggagcaacctgctctatfacagaccaactacctgggtggtggtgccatgatgattccattgtgggtttctgagt
cccttcaacatgatcctgggaggaaatcgtggtggtggtgttcacagggtttgtggtggcagccacaataaagacgtcttcgcccggatgaagaagcgtacccacgacgttcgttatggtggtcatgttggcgagctatttctatctccatgtttggaggatcgtgtttgtgtttg
gcattactttcttctgtgtgatgtttatccatgcgttgagacttcggaacctcaagaacaaactggagaataaaatggaaggaataggtttgaagaggacaccgatgggcattgtctggtggtccctagaacagcaggagaaggcatcaacagactcactgactatcagcaaat
gaaggaataa