



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

COMMD10-COMM domain containing 10 Gene

Reference: PTXBC005179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTD002; COMM domain-containing protein 10; COMM domain containing 10		
For research use only.			

GENE DESCRIPTION	COMM domain containing 10
GENE ID	51397
ORIGIN SPECIES	Human
ACCESSIONS	BC005179
VECTOR	pDONR223
NCBI SYMBOL	COMMD10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggtccccgcggcgctgactctcggagagccccagcatgagaagaacagctgctactgataaatgcaatagatcacaggaaatttccacggtgtctactcggattctcaaaaattccacctgaaggctgagagcagtttcagtgaagaagaggaaaaaaattcaagccgca-
ttttctctagagaanaacagatcttccactagtcttctgaacaataatcatatttttagaacagcagctgtatcacaaattgaagccagcagcttgcagcagcaattagagaacattcatctagacaagaacaagctgaagcattgtgtcaatcgtgtcttctatgggtcaagaacagctgaaag
tcccgccagagaaattctgctcctgttaagctagagaccgttgatggcagcttcaacttcagatggctcactcgtcctgaagcaaaaactaaattcctcaagctgtgttacaactcggagtgaacaatgaagattcaagacctgaggagaagttctgtggaattctgcacaaaggagttgttt
gggttcataaagacataacaagcagcagctgtgattcctctataga