

UTP11L-UTP11-like, U3 small nucleolar ribonucleoprotein, (yeast) Gene

Reference: PTXBC005182

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UTP11L; CGI-94; CGI94; U3 snoRNA-associated protein 11; UTP11-like protein; UTP11-like, U3 small nucleolar ribonucleoprotein; UTP11, small subunit processome component homolog (S. cerevisiae)		

For research use only.

GENE DESCRIPTION	UTP11-like, U3 small nucleolar ribonucleoprotein, (yeast)
GENE ID	51118
ORIGIN SPECIES	Human
ACCESSIONS	BC005182
VECTOR	pDONR223
NCBI SYMBOL	UTP11L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcggttttcggaaggcggcgttaagtcggcgcgggaacacagagagcgaagccagcctgcttgcgaaacatctgggctgctggagaaaaagaaagattacaaactctgcagatgactaccgtaaaaaaagaataacctcaagctcttcggaagaaggctctga
aaaaaaccagatgaattctactacaaatgactcgggttaactcaggatggagtacatattttaaggagactaaggagaagtaacccagaacaactaaagctgatgagaactcaggagctcaaatatataagaatgaagagggttcagaagctaagaaatcgaagactaaat
cagagctccatctgctggatttcaggggaagcaacagaacagcatgtgtttttgacacaaaaaggagttgaacagtttgatgtcgcaactcacctgcaaacagccccggagctagtcgacagagctcttaataaggccaggatagagaccttcgagaagaaaaagtgaaagg
agttaccaatcagactggacttaagcggatagctaaagaaggcaaaagcagataactgcctgacacagcggattgaacgagagaagaattgttcgttattgtcagaaattcaaacacgcaagatcttatgataaaactcagaagtgaaggtgaagaagaacggtgaactcc
ccagctattataaatttcagagtcgtcgaaacgttga