

UFM1-ubiquitin-fold modifier 1 Gene

Reference: PTXBC005193

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BM-002; C13orf20; ubiquitin-fold modifier 1		

For research use only.

GENE DESCRIPTION	ubiquitin-fold modifier 1
GENE ID	51569
ORIGIN SPECIES	Human
ACCESSIONS	BC005193
VECTOR	pDONR223
NCBI SYMBOL	UFM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgaagggtttccttaagatcacgtgacgtcgaccacggctgccgtacaaagtactcagtggtctgaaagtacaccttcacagcagctctaaagttgcagcagaagaattaaagttcctgctgcaacaagtgaattattaccaatgatggaataggaataaatcctgcacagactgctggaaatgttttctaaaacatggttcagaactcggattattcctagagatcgtgtggaagtgttaa