

MYOZ2-myozenin 2 Gene

Reference: PTXBC005195

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C4orf5; CMH16; CS-1; myozenin-2; FATZ-related protein 2; calcineurin-binding protein calsarcin-1; muscle-specific protein; myozenin 2		

For research use only.

GENE DESCRIPTION	myozenin 2
GENE ID	51778
ORIGIN SPECIES	Human
ACCESSIONS	BC005195
VECTOR	pDONR223
NCBI SYMBOL	MYOZ2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctatcacataatactatgatgaagcagagaaacagcaagcaacccatcatgaaggaagtcctggaatgatgttgatggcatggacgtgggcaaaaagtcagcatccccagagacatcatgttgaagaattatccatctcagtaaccgtggtgccaggctatttaagatgcg
tcaaagaagatctgacaaatacacatttgaaaattccagatcaatctagagcacaataaatcacagattgctatgcagaatgggaaagtggatggaagtaacttgaaggtggtcgcagcaagcccccttgactctctcccaacacccagatccacgaagccctccaaatccagacaa
cattgctccaggatattctggaccactgaaggaaattcctctgaaaaattcaacaccacagctgtccctaaagtaactatcaatctccctgggaacaagccattagcaatgatccggagcttttagaggctttatatctaaactttcaagcctgaaggaaaggcagaactgctgattacaggag
ctttaacaggggtgccacacatttgagggtttgaaaaagcatcaagaatggttaattaaagtccagatttgagctactattgtaacagatccaggtttatgtcctttgtcaatcccccttctggcagacggctcttaataggactcctaaggatggatctgagaatattcctatgtgat
aacaaccgaacctacagatgataccactgtaccagaatcagaagacctatga