



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KLK2-kallikrein-related peptidase 2 Gene

Reference: PTXBC005196

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KLK2A2; hGK-1; hK2; kallikrein-2; glandular kallikrein 2; glandular kallikrein-1; kallikrein 2, prostatic; tissue kallikrein-2; kallikrein related peptidase 2		

For research use only.

GENE DESCRIPTION	kallikrein-related peptidase 2
GENE ID	3817
ORIGIN SPECIES	Human
ACCESSIONS	BC005196
VECTOR	pDONR223
NCBI SYMBOL	KLK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggagacctgtcttctcatcgcccttctctgtggggctcactgggtgcgctgccctcatccactctgcgaattctgggagagcttggagatgtgagaagaactcccaaccttggcaggctgctgtgtacagctcatggatgggcacactgtgggggtgtcctgttgcacccccagtggtgtct
cacagctgcccaattgctaagaagaatagcaggcttggctgggtggcacaacactgtttgagcttgaagacacagccagagggtccctgttcagccacagcttccacacccgtctcataatatgagctcttgaagcatcaaaagctttagaccagatgaagactccgacatgacct
catgctgtcttgcgtctcagagccttccaaagatcacagatgttttgaaggtctctggcctgccacccagagacagcatctggggacacactgtcagcctcagcgtctgggcagatcgaaccagaggaagtcttggccccagaggtcttcagttgtgtgagccctcatctcttccaal
gacatgtgtgtctagagcttactctgagaaggctgacagagctcatgttctgtctggcctgtgcagagctgtgaagacacttctggggctgattctggggccactgtctgtaattgtgtgtcttcaggatcatcatcatctggggcccttgagccattgtccctgtcctgaaaagcctgtgtgta
caccaaggctgtgcttaccggaaagtgtgatcaagacacacatctcagccacccccctga