

ORMDL1-ORM1-like 1 (*S. cerevisiae*) Gene

Reference: PTXBC005200

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ORM1-like protein 1; adoplín-1; ORMDL sphingolipid biosynthesis regulator 1		
For research use only.			

GENE DESCRIPTION	ORM1-like 1 (<i>S. cerevisiae</i>)
GENE ID	94101
ORIGIN SPECIES	Human
ACCESSIONS	BC005200
VECTOR	pDONR223
NCBI SYMBOL	ORMDL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacgttggagttgccacagtgagtgaaatccaaataccgtgtcatgaacagccgggtatgtgctgacatagcattgggagttggcttcatattgtcttactcagcattcccttctcagtggttctgttggactttaacaaatattacataatctgggatgtacgtatttg
catgcagtgaagggaacaccttccgaaactcctgaccaggttaagcaaggctcctaactcattgggaacaactggactatggagtagcatttacatcttcacggaagttttcacaatttcccaataattctatatttctggcaggtttctatcgaagtatgatccaactcacttcactcctaaca
cagcttctctcctgagtgactaattcccaaatgccacaactacatgggttcggatctttggaattaataagtattga