

## NIPSNAP3B-nipsnap homolog 3B (C. elegans) Gene

Reference: PTXBC005202

|                               |                                                                        |                     |                |
|-------------------------------|------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | FP944; NIPSNAP3; SNAP1; protein NipSnap homolog 3B; nipsnap homolog 3B |                     |                |
| <b>For research use only.</b> |                                                                        |                     |                |

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | nipsnap homolog 3B (C. elegans)                                                                   |
| <b>GENE ID</b>          | 55335                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC005202                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | NIPSNAP3B                                                                                         |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtctgttctcagaagcggcctgaccaagcgcttgccctacggagctcgccctcagggtgtgttcattctttgctacggccctagacaatcagatggaacgttctatgaatttcgtacttattaccttaaacctcaaatatgaatgcgttcaggaataatctaagaaaaacattcatctcgg  
acctcttactctgaattggttgattctggagtgtagaatttggaggcagaacgaataaagtgttcatatttggagtgatgataatttgcctcagctgagctgaagttcggaaagccttagccaactgtaaggaaatggcaagaacaatctatctccaaatttggtcgcattgataaacaagagac  
ggaaattacttacctgataccatggtccaattagaaaaagcctccaaaagaggagctctgaactagctgttttcagatgaacctgggtggccagctctgtgggtgatgcattgaaagagcaatgaatgccatgtcaatttaggctacacaaaagtagttggtgtttccacacagaat  
ggagaactcaacagagttcatgttcttgggtgaatgagagtgagatagctgtcagctgggagacataagtcctcagagatccagagttgtggcggctgttcgggaagtgfcaactacgtgttctcagcagaatgtcttctgattcctgcatttccattgaaatag