

FXYD3-FXYD domain containing ion transport regulator 3 Gene

Reference: PTXBC005238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAT8; PLML; FXYD domain-containing ion transport regulator 3; chloride conductance inducer protein Mat-8; mammary tumor 8 kDa protein; phospholemman-like protein; FXYD domain containing ion transport regulator 3		

For research use only.

GENE DESCRIPTION	FXYD domain containing ion transport regulator 3
GENE ID	5349
ORIGIN SPECIES	Human
ACCESSIONS	BC005238
VECTOR	pDONR223
NCBI SYMBOL	FXYD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagaaggtgacctggcctgtgtgtcctggcaggcttctgtcctggagccaatgacctaagataaaaacagtccttctactatgactggcacagcctccaggttgccgggctcatctgcgctggggttctgtgcgccatgggcatcatcatcgtcatgagtgcaaaatgca
aafgcaagtttggccagaagtccggtcaccatccaggggagactccacctctcatcaccacaggtcagcccaagctga