

HAX1-HCLS1 associated protein X-1 Gene

Reference: PTXBC005240

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HCLSBP1; HS1BP1; SCN3; HCLS1-associated protein X-1; HAX-1; HCLS1 (and PKD2) associated protein; HS1 binding protein; HS1-associating protein X-1; HS1-binding protein 1; HSP1BP-1; HCLS1 associated protein X-1		

For research use only.

GENE DESCRIPTION	HCLS1 associated protein X-1
GENE ID	10456
ORIGIN SPECIES	Human
ACCESSIONS	BC005240
VECTOR	pDONR223
NCBI SYMBOL	HAX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctcttgcctcctcggggcttttgccttctgacctcggagccacagagatcccttttggaggggactcagatgaagatgatgatgaggaagaagaagaagggggctcatggggccgtgggaaccaaggttcctagccctcagcaccctcaggggaatt
tggcttcggcttcagcttcagccagaggaggaggtacgttccacgataactcggcttgcctgacacgtacgagattcaatagcatcttcagcgatgaggggctggaccttgcctccatcctcctgaacttccaggtcctgagtcagagacacctggtagagactacggaggga
cagacacttcgggactcaatgcttaagtatccagatgacccagccaggaatcttggggggcttggagagtgatgcaagaagtgaatcccccaaccagcaccagactggggctccagaggccattcataggttgcctatggaccccatcctagaaccagag
aggacaatgatcttgattccaggttccagagggtcttggccgggtctacagccagcccaatcctattcaagagcatctctgtgaccaagatcactaaaccagatgggatagtgaggagcgcgggactgtgtgacagtggggccggacagactacagtaaccgac
acgaagcagatagcagtcctaggggtgatccagaatcaccaagacctccagccctggatgacgttccatcctgcacttattcctgggacgttgggtccgggtcccgtag