

TSPAN8-tetraspanin 8 Gene

Reference: PTXBC005246

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CO-029; TM4SF3; tetraspanin-8; transmembrane 4 superfamily member 3; tspan-8; tumor-associated antigen CO-029; tetraspanin 8		

For research use only.

GENE DESCRIPTION	tetraspanin 8
GENE ID	7103
ORIGIN SPECIES	Human
ACCESSIONS	BC005246
VECTOR	pDONR223
NCBI SYMBOL	TSPAN8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggtgtgagtgccgtataaaatattctatgtttacctcaactctgttctggtctatgtgtatcttgatcctagcattagcaaatgggtacgaataagcaatgactctcaagcaattttggtctgaagatgtaggctctagctcctacgttgctgtggacatatgattgctgtaggcgcca
tcacatgattctgggcttccctggcatgctgctgtctataaaagaagtcgctgcatgcttctgttttttcatagccttgcctctgcatgctcctgcaggtggcgacaggtatcctaggagctgttttcaaatcaagtctgacgcattgtgaatgaaactctctatgaaaacacaaagcttttg
agcgccacaggggaagtgaaaaacaattccagggaagccataatgtgttcaagaagagtttaaatgctgctgtttgttcaatgagctgctgattggggaataatttcaacactatcctgaattatgtgcctgtctagataagcagagaccatgccaaagctataatggaaaacaagtta
caaagagacctgtattcttataaaagactcttgcaaaaaatttgattatagtattggaaatagcattggactggcagttattgagatactgggttgggttttctatgtgctctgtattgccagatcggaacaaatga