

COX7A2L-cytochrome c oxidase subunit VIIa polypeptide 2 like Gene

Reference: PTXBC005251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	COX7AR; COX7RP; EB1; SIG81; cytochrome c oxidase subunit 7A-related protein, mitochondrial; COX7a-related protein; cytochrome c oxidase subunit VII-related protein; cytochrome c oxidase subunit VIIa polypeptide 2 like; cytochrome c oxidase subunit VIIa-re		

For research use only.

GENE DESCRIPTION	cytochrome c oxidase subunit VIIa polypeptide 2 like
GENE ID	9167
ORIGIN SPECIES	Human
ACCESSIONS	BC005251
VECTOR	pDONR223
NCBI SYMBOL	COX7A2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtactacaagtttagtggcttcacgcagaagttggcagagcatgggcttcggagcctatagcccgacgggattaaagcctgtgtttccacagaagcaccacatcatattgccacacaaactaaactgacctccgattccacagtgtatgattatgctgggaaaaacaaagttccag
agctacaaaagttttccagaaaagctgatgtgtgtcccgtctacctgaaacgagggcctgacctgaacaaatgctttaccggaccacatggcgctgactgtggaggggaccatctactgctgatgcctctcatatggcttcgcagcccaaaaacaatga