

KLHL20-kelch-like 20 (Drosophila) Gene

Reference: PTXBC005253

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KHLHX; KLEIP; KLHLX; kelch-like protein 20; Kelch motif containing protein; kelch-like ECT2-interacting protein; kelch-like protein X; kelch like family member 20		

For research use only.

GENE DESCRIPTION	kelch-like 20 (Drosophila)
GENE ID	27252
ORIGIN SPECIES	Human
ACCESSIONS	BC005253
VECTOR	pDONR223
NCBI SYMBOL	KLHL20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggaaagccaatgccaggtgtaccaacattcgaccaggagagactggaatggatgaacaagccgtgcaccctggagacccaacaactgccagaaggggtcccaacctgcccgcattatctcagacaagcaccctcgacaaccttgaagtattaaacct
tctgagaagcaccgggagctatgtgatgtgtgctagtgtggcgccaagaagatatatgccatcgagtcattttgtcagcctgtatgccctactccgagctatgtttacaggagaattggcagagagccgtcagacagaagtagtgcagacattgacgagagggctatggaatt
actgattgactttgctatacctccagataacagtagaagagggaatgttcagactcttctgccagctgcttgcctccagctggcagaatacaggagcctgctgtgaattctaaagagacaattagatccttctaactgctgggcatcgggcttttctgacacacattcatgtcgt
gagttgctaaggatagcagacaagttcacccaacatacttcaagaggtgagttgcactgtgttccaatgcacccattgttgagagcaatatgggtcacattctgtaataaacactgcctgggtatttctctgctga