

CLEC2B-C-type lectin domain family 2, member B Gene

Reference: PTXBC005254

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AICL; CLECSF2; HP10085; IFNRG1; C-type lectin domain family 2 member B; C-type (calcium dependent, carbohydrate-recognition domain) lectin, superfamily member 2 (activation-induced); C-type lectin superfamily member 2; IFN-alpha-2b-inducing-related protei		

For research use only.

GENE DESCRIPTION	C-type lectin domain family 2, member B
GENE ID	9976
ORIGIN SPECIES	Human
ACCESSIONS	BC005254
VECTOR	pDONR223
NCBI SYMBOL	CLEC2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgaccaacataaaaagtgtttataattgttggtgttttaataacaataattattactctgatagttaaactaactcgagattctcagagtttatgccctatgattggattggttccaaaacaaatgctattatttctaaagaagaaggagattggaattcaagtaatacaactgttccact
caacatgccgacctactataattgacacatagaagaacgaattttcttaggcgtataaatgcagttctgatcactggattggactgaagatggcaaaaatcgacaggacaaatgggtagatggagctacattaccaaatcgttggcatgagaggagtgaggatgtgcctactca
gcgatgatggtgcagcaacagctagatttacaccgaaagaaatggatttcagggaaagaatacactaa