

MSMB-microseminoprotein, beta- Gene

Reference: PTXBC005257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HPC13; IGBF; MSP; MSPB; PN44; PRPS; PSP; PSP-94; PSP57; PSP94; beta-microseminoprotein; immunoglobulin binding factor; prostate secreted seminal plasma protein; prostate secretory protein of 94 amino acids; seminal plasma beta-inhibin; microseminoprotein		

For research use only.

GENE DESCRIPTION	microseminoprotein, beta-
GENE ID	4477
ORIGIN SPECIES	Human
ACCESSIONS	BC005257
VECTOR	pDONR223
NCBI SYMBOL	MSMB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgttctcctggcagcgtgtgatctttgccacctcgtgactttatgcaatgcatcgtctatttcatacctaaggggagttccaggagattcaaccaggaaatgcatggtctcaaaggaaacaaacaccaataaaactggagtggcagactgacaactgtgagacatgcactgtacgaacagaaatttcattgtgacccctgtttctacacctgtgggttatgacaaagacaactgccaaagaatctcaagaaggaggactgcaagtatatctgtgtggagaagaaggacccaaaaagacctgttctgtcagtgatggataatctaa