

INPP4B-inositol polyphosphate-4-phosphatase, type II, 105kDa Gene

Reference: PTXBC005273

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	type II inositol 3,4-bisphosphate 4-phosphatase; inositol polyphosphate 4-phosphatase II; 4-phosphatase II; inositol polyphosphate-4-phosphatase, type II, 105kDa; inositol polyphosphate-4-phosphatase type II B		

For research use only.

GENE DESCRIPTION	inositol polyphosphate-4-phosphatase, type II, 105kDa
GENE ID	8821
ORIGIN SPECIES	Human
ACCESSIONS	BC005273
VECTOR	pDONR223
NCBI SYMBOL	INPP4B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaattaagaggaggcatcagaagaaggcagcactttcttctacagccaggccaatgatccggggactgtcagttcacaagatccagaagactccaaatgaaccgcagttggaattcatccttgattgaatcaagaattaagagattag