

APOA2-apolipoprotein A-II Gene

Reference: PTXBC005282

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Apo-AII; ApoA-II; apoAII; apolipoprotein A-II; apolipoprotein A2		

For research use only.

GENE DESCRIPTION	apolipoprotein A-II
GENE ID	336
ORIGIN SPECIES	Human
ACCESSIONS	BC005282
VECTOR	pDONR223
NCBI SYMBOL	APOA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgtcgcagcaactgtgctactctcaccatctgcagcctgaaggagctttggttcggagacaggcaaggagccatgtgtggagagcctggtttctcagtacttcagaccgtgactgactatggcaaggacctgatggagaaggtcaagagccagagctcagccgag
gccaaagtctactttgaaaagtcaggagcagctgacacccctgatcaagaaggctggaacggaactggttaactcttgagctattctgtggaacttggaacacagcctgccaccagtgga