

PRM2-protamine 2 Gene

Reference: PTXBC005303

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT94.2; protamine-2; cancer/testis antigen family 94, member 2; sperm histone P2; sperm protamine P2; testicular secretory protein Li 40; protamine 2		

For research use only.

GENE DESCRIPTION	protamine 2
GENE ID	5620
ORIGIN SPECIES	Human
ACCESSIONS	BC005303
VECTOR	pDONR223
NCBI SYMBOL	PRM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtccgataccgctgaggagcctgagcgaacgctgcacgaggtgtacaggcagctgtcatgggcaagagcaaggacaccacggccaagaggagcaagggtgagccgggagcacgtcgaggtctacgagaggacccatggccagtctactataggcgacactgc
tctcgaaggaggctgcaccgatccacaggcgagcagcatcgtctctcagaaggcgcaaaagacgtctctgcaggcaccggaggagggatcgagaggtcgagaaccaggaagagaacatgcagaaggcactaa