

GLRX-glutaredoxin (thioltransferase) Gene

Reference: PTXBC005304

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GRX; GRX1; glutaredoxin-1; TTase-1; thioltransferase; thioltransferase-1; glutaredoxin		

For research use only.

GENE DESCRIPTION	glutaredoxin (thioltransferase)
GENE ID	2745
ORIGIN SPECIES	Human
ACCESSIONS	BC005304
VECTOR	pDONR223
NCBI SYMBOL	GLRX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcaagagtttgtgaactgcaaaatccagctgggaaggtgtgtgtgtcatcaagccacctgccgtactgcaggaggcccaagagatcctcagtcattgccatcaacaaggctctggaattgtcgatatacagccaccaaccacactaacgagattcaagattattgc
aacagctcacgggagcaagaacgggtgcctcgagctcttatttgtaaaagtgtatagcggatgcagtgatctagtctcttgcacagagtggggaactgctgacgcggctaaagcagattggagctctgcagtaa