

COL4A6-collagen, type IV, alpha 6 Gene

Reference: PTXBC005305

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXDELq22.3; DELXq22.3; DFNX6; collagen alpha-6(IV) chain; collagen IV, alpha-6 polypeptide; collagen of basement membrane, alpha-6; collagen, type IV, alpha 6; dJ889N15.4 (Collagen Alpha 6(IV)); collagen type IV alpha 6 chain		

For research use only.

GENE DESCRIPTION	collagen, type IV, alpha 6
GENE ID	1288
ORIGIN SPECIES	Human
ACCESSIONS	BC005305
VECTOR	pDONR223
NCBI SYMBOL	COL4A6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttataacaagttgtgctgctcctgtttagctgtgacctgacgaggaactggcagcagcgggagagaagtcttatgaaagccatgtggggccaggactgcagtgggagctgtcagtgcttttctgagaaaggagcgagacacaacttcagctactgaatgatatggctggccggctctaccatttcagtgagggttcgcctaatactcttttga