



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KLK3-kallikrein-related peptidase 3 Gene

Reference: PTXBC005307

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APS; KLK2A1; PSA; hK3; prostate-specific antigen; P-30 antigen; gamma-seminoprotein; kallikrein-3; semenogelase; seminin; kallikrein related peptidase 3		

For research use only.

GENE DESCRIPTION	kallikrein-related peptidase 3
GENE ID	354
ORIGIN SPECIES	Human
ACCESSIONS	BC005307
VECTOR	pDONR223
NCBI SYMBOL	KLK3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggctccggctgtcttctccacctgtctcgtgacgttgattggctgcgcgccctcatctgtctcgallttggggagactgggaagtgcgaagaacatcccaacctggcaggctgtgtgtgcctctgtggcagggcagctgcggcggtgttctgtgtgcacccacagtggtctctc
acagctgcccaactgcatcagggaacaaaaggctgattctgtcgttgctggccacagctgtttctatctggaagacacaggccaggtatttcaggctgacccacagcttccacacccgcctctacgatattgagctcctgaagaatcattctcaggccaggtgtgtagtactccagccacagctcat
gtctctcgcctgtcagagccttgcagagctcagggtgctgtgaagctcatggacctgcgccacaggagacagacatctggggacacactgtctagcctcaggctggggcagcatgaaccagagagattcttgaccacaaagaactcaagtgtgtgacctccatgtatttccaatga
cgtgtgtgcgaagtccacctcagaaggtgaccaagttcatctgtgtgtgacctgcggacgtcgagacgtggggccaaaggcgaagctgtgggtgattctggggcccatctgtgtgattgtgtgtcttaagttatcatcgtatcgttcgggcagtgtaacctgtgacctgccgaaggccttccctgt
acacaaagtgctgcttaccgggaattgcgacaaagacacatcattgcgcccaacctcga

More info: Please contact us for any question: contact@proteogenix.fr