



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TTR-transthyretin Gene

Reference: PTXBC005310

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CTS; CTS1; HEL111; HsT2651; PALB; TBPA; ATTR; carpal tunnel syndrome 1; epididymis luminal protein 111; prealbumin, amyloidosis type I; thyroxine-binding prealbumin

For research use only.

GENE DESCRIPTION	transthyretin
GENE ID	7276
ORIGIN SPECIES	Human
ACCESSIONS	BC005310
VECTOR	pDONR223
NCBI SYMBOL	TTR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcctcaatgctgctgctcctctgctctgctgagctgattgttgctgagggctggccctacggggcaccggctgaatccaagtgtctctgatgtgtcaaaagtctgatgtctgccgaggcgactctgccatcaatgtggcgctcatgtgttcagaaaaggctgctgatgacacctgggagcctgtgctcctctggaaagacacagctgagctgctgagagctgctgagctctcaacactgaggaagattgtgagaaggatatacaaaatggcacaatagacaccaaaatcttactggaaggcgctggcatctcccaatcctatgacatgctcagagagtgctgtattcacagccaacgactcgggccccccctctgacacattgcccgcctctgtagcctcattccaccacggcgtctgctcaccacatcccaaggaaatga