

GYPA-glycophorin A (MNS blood group) Gene

Reference: PTXBC005319

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD235a; GPA; GP _E erik; GPSAT; HGpMiV; HGpMiXI; HGpSta(C); MNS; PAS-2; glycophorin-A; MN sialoglycoprotein; Mi.V glycoprotein (24 AA); erythroid-lineage-specific membrane sialoglycoprotein; glycophorin A (MN blood group); glycophorin A, GPA; glycophorin Erik		

For research use only.

GENE DESCRIPTION	glycophorin A (MNS blood group)
GENE ID	2993
ORIGIN SPECIES	Human
ACCESSIONS	BC005319
VECTOR	pDONR223
NCBI SYMBOL	GYPA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatggaaaaataatcttcttattactattgtcagcaattgtgagcatatcagcatttaagtaccactgaggtggcaatgcacacttcaacctctcttcagtcacaaagagttacatctcatcacagacaaatgatacgcacaaacgggacacatatgcagccactcctagagctcatgaagttc
agaaattctgttagaactgtttaccctccagaagaggaaacccggagaagggtacaactgcccatttctctgaaccagagataacactcattatctttgggtgatggctggtgttattggaacgatctcttaatttctacggtattcgcgactgataaagaaaagcccatctgatgtaa
acctctccctcactgacacagacgtgcctttaagttctgttgaatagaaaatccagagacaagtgaatga