

TNNC2-troponin C type 2 (fast) Gene

Reference: PTXBC005323

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	troponin C, skeletal muscle; troponin C type 2 (fast); troponin C2, fast; troponin C2, fast skeletal type		

For research use only.

GENE DESCRIPTION	troponin C type 2 (fast)
GENE ID	7125
ORIGIN SPECIES	Human
ACCESSIONS	BC005323
VECTOR	pDONR223
NCBI SYMBOL	TNNC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggaccagcaggtgagggcaggtctacacagcgaagagatgctgctgagttcaaggctgctttgacatgtttgatgctgatggtggtgggacacagcgtcaaggagttgggcacgggtgatgaggatgctgggacagacacccaccaaggaggagctggacgccatca
tcgaggaggtggatgaggacggcagcggcaccatcgacttcgaggagttcttggatcatggtgctgagatgaaagaggacgcgaaagggaagagcaggaggagctggccgagtgcttccgcatcttcacaggaatgcagacggctacatcgaccggaggagctggctga
gatttcaggccctccgggagcagtgactgacgaggagatcgaatctctgatgaaagacggcgacaagaacaacgacggccgattgacttcgacgagttctgaagatgatggaggcgtgcagtaa