

RNASE1-ribonuclease, RNase A family, 1 (pancreatic) Gene

Reference: PTXBC005324

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAC1; RIB1; RNS1; ribonuclease pancreatic; HP-RNase; RIB-1; RNase 1; RNase A; RNase upI-1; ribonuclease 1; ribonuclease A C1; ribonuclease, RNase A family, 1 (pancreatic); ribonuclease A family member 1, pancreatic		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, 1 (pancreatic)
GENE ID	6035
ORIGIN SPECIES	Human
ACCESSIONS	BC005324
VECTOR	pDONR223
NCBI SYMBOL	RNASE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctggagaagtctcttgcggctccttctgttgcctgatactgtggtgctgggtgggtccagccttccctgggcaaggaaatccgggcaagaattccagcggcagcatatggactcagacagttccccagcagcagctccactactgtaacaaatgatgagcgccg
gaatatgacacagggcggtgcaaaccaagtgaacaccttgtgcacgagccctggtagatgtccagaatgtctgtttccaggaaaaggtcacctgcaagaacgggcaggcgcaactgctacaagagcaactccagcatgcacatcacagactgccgcctgacaaacggctccaggtac
cccaactgtcataccggaccagcccgaaggagagacatcatgtggcctgtgaaggaggcccatatgtgccagtcactttgatgcttctgtggaggactctaccta