

RPL27A-ribosomal protein L27a Gene

Reference: PTXBC005326

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L27A; 60S ribosomal protein L27a; ribosomal protein L27a		

For research use only.

GENE DESCRIPTION	ribosomal protein L27a
GENE ID	6157
ORIGIN SPECIES	Human
ACCESSIONS	BC005326
VECTOR	pDONR223
NCBI SYMBOL	RPL27A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccatccagactgaggaagaccggaacttagggccacgtgagccacggcgcgcagtaggcaagcaccggaagcaccgcccggcgccggtatgctgggtctgcatcaccacggatcaactcgacaaataccaccaggctacttgggaaagttggtatgaa
gcattaccacttaagagggaaccagagcttctgccaactgtcaaccttgacaaattgtggactttggtcagtgacagacacgggtgaatgctgctaaaaacaagactgggctgctcccatcattgatgtggtcgatcgggctactacaaagtctgggaaagggaagctcccaagc
agcctgtcatcgtgaaggccaaattctcagcagaagagctgaggagaagattaagagtgttggggggcctgtgtcctgtgcttga