

CETN2-centrin, EF-hand protein, 2 Gene

Reference: PTXBC005334

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CALT; CEN2; centrin-2; caltractin (20kD calcium-binding protein); centrin, EF-hand protein, 2; centrin 2		

For research use only.

GENE DESCRIPTION	centrin, EF-hand protein, 2
GENE ID	1069
ORIGIN SPECIES	Human
ACCESSIONS	BC005334
VECTOR	pDONR223
NCBI SYMBOL	CETN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctccaactttaagaaggcaaacatggcatcaagttctcagcgaaaaagaatgagccctaagcctgagccttactgaagagcaaaagcaggagatccgggaagcttttgatcttttcgatgcggatggcaactggcaccatagatgttaagaactgaaggtggcaatgagggccctgg
gctttgaacccaagaaagaagaataaagaaatgataagtgaattgataagggaaggaacaggaataatgaactttggtgactttttaactgtgatgacccagaaaatgtctgagaagatactaaagaagaatcctgaagcttcaagctttgatgatgaaactgggaagatttcg
ttcaaaaatctgaaacgcgtggccaaggagttgggtgagaacctgactgatgaggagctgcaggaaatgattgatgaagctgacgagatggagatggagaggtcagtgagcaagagttcctgcgcatcatgaaaaagaccagcctctattaa