

PSMA4-proteasome (prosome, macropain) subunit, alpha type, 4 Gene

Reference: PTXBC005361

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HC9; HsT17706; PSC9; proteasome subunit alpha type-4; macropain subunit C9; multicatalytic endopeptidase complex subunit C9; proteasome (prosome, macropain) subunit, alpha type, 4; proteasome component C9; proteasome subunit HC9; proteasome subunit L; pro		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) subunit, alpha type, 4
GENE ID	5685
ORIGIN SPECIES	Human
ACCESSIONS	BC005361
VECTOR	pDONR223
NCBI SYMBOL	PSMA4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcgaagatagactccaggaccactatatttccagaaggtcgttataccaagtgaatatgcatggaagctattggacatgcaggcacctgtttgggaatttagcaaatgatggtgtttgcttgcagcagagagacgcaacatccacaagcttctgatgaagtctttttctgaaa
aaatttataaactcaatgaggacatggcttgcagtggtgcaggcagataactctgatgctaattgttctgactaatgaactaaggctcattgctcaaaaggattattacagatcaggagccaataaccttgcagcaggttggttacagcgtgtgtgatatacaaacagcttatacacaatttggaggaa
aacgtccctttgtgtttcattgctgtacattggctgggataagcactatgctttcagctctatcagagtgaccctagtggaattacggggatggaaggccacatgcattggaaataatagcgtgcagctgtgtcaatgttgaacaagactataaagaaggagaatgacctgaagtca
gcacttgccttagctatcaagactaataagaccatggatgttagtaactctctgctgaaaaagtgaattgcaactaacaagagagaatggaagacagtaacagagttctcaacaaaaagtgagcagttgatcaaaaacacgaggagaagaagccaaagctgagc
gtgagaagaagaaaaaacagagaagaaaaagataaatag