

MRPS18B-mitochondrial ribosomal protein S18B Gene

Reference: PTXBC005373

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf14; HSPC183; HumanS18a; MRP-S18-2; MRPS18-2; PTD017; S18amt; 28S ribosomal protein S18b, mitochondrial; 28S ribosomal protein S18-2, mitochondrial; MRP-S18-b; S18mt-b; mitochondrial ribosomal protein S18-2; mrps18-b; mitochondrial ribosomal protein S		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S18B
GENE ID	28973
ORIGIN SPECIES	Human
ACCESSIONS	BC005373
VECTOR	pDONR223
NCBI SYMBOL	MRPS18B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgctgtatataaacaccgtgctgaggcggcttctatgtctatctcttccagggttcacagagttcaggttccctccagactcttgcacaaagctccctctgagggaagattcttctctcagttccatttctcttataagatgagccctggaatatctggaatcagaagaa
taccaggagcgataggttctcgcccgctctggcctgactaccgcgaaccacaagggtgtgtacccccacagcggactcggagacatgtattctcgggaataaagtgttgggaatccctgccccatctctcgagatcacaggttgcattttaggaacgtgaagctcttgag
caatttctcgcgccacacgggtatcttctatgtctccatacacaggagttgtgtgaagcagcacaagcgttgacccagccatccagaagccaggatcatgtctctcatttaccacatccccaggttgaaccacgggacctgacttacttctcatgggctgtgagtg
ctactccgcagccccaccctgtctcaagtaccctgtacccatggtacaactgaaacagccaccggagagagaactgtctgccttcgcccgtttaccagggtcatctccaagaagagagtgccccccacctgagtcattgccaagatgccccctagaacaccagcggaa
gcctctccactgggcagacaggccctcagagtgctctgtag