

APOA1-apolipoprotein A-I Gene

Reference: PTXBC005380

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	apo(a); apolipoprotein A-I; apo-AI; apolipoprotein A1		
For research use only.			

GENE DESCRIPTION	apolipoprotein A-I
GENE ID	335
ORIGIN SPECIES	Human
ACCESSIONS	BC005380
VECTOR	pDONR223
NCBI SYMBOL	APOA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagctgcggtgctgacctggccgtgctctctgacggggagccagctggcattctgagcagcaagatgaacccccagagccctgggacgagtgaggacctggccactgtgtacgtggatgtgctcaagacagcggcagagactatgtgccagtttgaaggctcc
gccttgggaaaacagctaaacctaaagctccttgacaactgggacagcgtgacctccacctcagcaagctgcgcgaacagctcgccctgtgacccaggagttctgggataacctggaaaaggagacagagggcctgaggcaggagatgagcaaggatctggaggaggtgaaggc
caaggtgcagccctacctggacgactccagaagaagtggcaggaggagatggagctctaccgccagaaggtggagccgctgcgcgcagagctccaagaggcgcgccagaaagctgcacgagctgcaagagaagctgagccactggcgaggagatgcgcgacccgcgcg
cgccccatgtgacgcgctgcgcgacgcatctggccccctacagcgacgagctgcgccagcgttggccgcgcgcctgaggctctcaaggagaacggcgccagactggccgagttaccacgccaaggccaccgagcatctgagcagctcagcgagaaggccaagccgcg
gctcaggacctccgaaggcctgctgcccgtgctggagagctcaaggctcagctctcagagagtgagtagtaagaagctcaacaccagtgga