

CETN3-centrin, EF-hand protein, 3 (CDC31 homolog, yeast) Gene

Reference: PTXBC005383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDC31; CEN3; centrin-3; CDC31 yeast homolog; EF-hand superfamily member; centrin, EF-hand protein, 3 (CDC31 homolog, yeast); centrin 3		

For research use only.

GENE DESCRIPTION	centrin, EF-hand protein, 3 (CDC31 homolog, yeast)
GENE ID	1070
ORIGIN SPECIES	Human
ACCESSIONS	BC005383
VECTOR	pDONR223
NCBI SYMBOL	CETN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtttagctctgagaagtgcctttagtgacacaaacagaggagaaaaaagaagagaactgtctgaggaacagaaacagaatgaagatgctttgaactatttgatacagacaaagatgaagcaatagattatcatgaattaaaggtggcaatgagagccttggggttgatga
aaaaaagctgactgaagattcttaagattatgacagagaagccacaggggaaatcaccttgaagattttaaagaagttgtgacagactggatattggaaagagatcccatgaagaaatactcaaggcatttaaactatttgatgatgattcaggtaaaataagcttgaggaaattgc
gacgtgtgctagagaattgggtgaaacatgagtgatgaagaacttcgagctatgatagaagaattgacaaagatggtgatggagaataaaccaagaggaggttcattgctattatgactggtgacatttaa