

PLA2G1B-phospholipase A2, group IB (pancreas) Gene

Reference: PTXBC005386

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLA2; PLA2A; PPLA2; phospholipase A2; phosphatidylcholine 2-acylhydrolase 1B; phospholipase A2, group IB (pancreas); phospholipase A2 group IB		

For research use only.

GENE DESCRIPTION	phospholipase A2, group IB (pancreas)
GENE ID	5319
ORIGIN SPECIES	Human
ACCESSIONS	BC005386
VECTOR	pDONR223
NCBI SYMBOL	PLA2G1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaactcttgtgctagctgtgctgctcacagtggccgcccgcacagcggcatcagccctcggccgtgtgagcagttccgcaaaatgatcaagtgcgtgatccggggagtgaccctcttggatacaacaactacggctgctactgtggcttgggggctcaggcaccctcgtgatgaactggacaagcaaaaacaagagtgatga