



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C11orf51-chromosome 11 open reading frame 51 Gene

Reference: PTXBC005393

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf51; APC15; HSPC020; anaphase-promoting complex subunit 15; anaphase promoting complex subunit 15		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 51
GENE ID	25906
ORIGIN SPECIES	Human
ACCESSIONS	BC005393
VECTOR	pDONR223
NCBI SYMBOL	C11orf51
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccaattgttcccctcaactctccctcgtgtgactgagacctgtgtgttaatactgcatgacacctgtgtgtggaagacagacagctgcagcagcaggaacagcagcatcaggctcgtgtcccaaaagatcgcgcgagaagaacacaacctgtgttcatttggcaagccagcttcagagca
ctatgatgacgaggaagaagaggatgatgaagatgatgaggatagtgaaggaactcagaggatgatgaggatagtgcaggacatgcagcagatgaatgactacaatgagtcaccgcgatgatggagaggtcaatgagtggtgacatggaagcgaacgaacaggatcaggaccagtga
tcatctag