

ANAPC13-anaphase promoting complex subunit 13 Gene

Reference: PTXBC005398

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	APC13; SWM1; anaphase-promoting complex subunit 13; cyclosome subunit 13; anaphase promoting complex subunit 13		

For research use only.

GENE DESCRIPTION	anaphase promoting complex subunit 13
GENE ID	25847
ORIGIN SPECIES	Human
ACCESSIONS	BC005398
VECTOR	pDONR223
NCBI SYMBOL	ANAPC13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagtgggtcagagagatggaaggatctggatttgattgatgcttgccgagagaagacaagctgccttatgaggatgtcgcaataccactgaatgagcttcctgaacctgaacaagacaatggtgcaccacagaatctgtcaagaacaagaatgaagtgacagacttag
ccttacagtacctccatgagaatgttcccccatggaaactga