

UCRC-ubiquinol-cytochrome c reductase complex (7.2 kD) Gene

Reference: PTXBC005402

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UCRC; HSPC051; HSPC119; HSPC151; QCR9; UCCR7.2; cytochrome b-c1 complex subunit 9; cytochrome C1, nonheme 7kDa protein; cytochrome c1 non-heme 7 kDa protein; ubiquinol-cytochrome c reductase complex (7.2 kD); ubiquinol-cytochrome c reductase complex 7.2 k		

For research use only.

GENE DESCRIPTION	ubiquinol-cytochrome c reductase complex (7.2 kD)
GENE ID	29796
ORIGIN SPECIES	Human
ACCESSIONS	BC005402
VECTOR	pDONR223
NCBI SYMBOL	UCRC
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccgcgacgttgactcgaagtggtactccctgctgttccgcaggacctccaccttcgacctcaccatcatcgtggcgctcatgttcttcgagcgcgccttcgatcaaggcgcgacgctatctacgaccacatcaacgaggggaagctgtggaacacatcaagcacaagtatga
gaacaagtag