

KIAA0101-KIAA0101 Gene

Reference: PTXBC005832

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA0101; NS5ATP9; OEATC; OEATC-1; OEATC1; PAF; PAF15; p15(PAF); p15/PAF; p15PAF; PCNA-associated factor; HCV NS5A-transactivated protein 9; PCNA-associated factor of 15 kDa; hepatitis C virus NS5A-transactivated protein 9; overexpressed in anaplastic thy		

For research use only.

GENE DESCRIPTION	KIAA0101
GENE ID	9768
ORIGIN SPECIES	Human
ACCESSIONS	BC005832
VECTOR	pDONR223
NCBI SYMBOL	KIAA0101
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcgcactaaagcagacagtggtccaggcacttacagaaaagtggtgctgctcagccccagaaaggtgcttggttctccacctctgccactaatcgacatcagttcatcgaggaaagctgaaataaatatgcaggagggaaccccggttgcgtgcgccaactccaagtgcaaaaaggaaattggagaattctttaggttgctccctaagattctgaaaaagagaatcagattctgaagaggcagggaagcagtggtcttaggaaaagcaagagaaaagcatgtcctttgcaacctgatcacacaaatgatgaaaaagaatag